

Safety Lancets

For all your point-of-care blood sampling needs



**Safety,
simplicity
and
comfort**



Efficiency in action

The adjustable-depth safety lancet

When you're ready to draw a sample, the last thing you want is to be searching around for the right depth lancet. Our adjustable-depth safety lancet has three depth settings to minimize skus and increase efficiency.

- Innovative rotating dial allows for three easy depth adjustments: 1.3 mm, 1.8 mm, 2.3 mm
- Needle remains safely covered, preventing accidental sticks
- Pressure activated lancet is simple and effective
- Compares to the Accu-Chek® Safe-T-Pro Plus

Item No.	Description	Pkg.
MPHADJDEPTH	23 gauge, Adjustable-Depth Safety Lancet	1,000/cs
MPHADJDEPTHZ	23 gauge, Adjustable-Depth Safety Lancet	200/bx



Blood sampling **simplified**

Meet every patient's need and every caregiver's preference

You perform hundreds of blood draws every week. You want them to be quick, painless and free of accidental needle stick injuries.

Our entire line of lancets is designed to be safe and user-friendly. No matter what style you prefer, we have what you need. And now we offer a new adjustable depth lancet to provide individualized patient care.



Push button
MPHSFTY28

Compares to
Accu-Chek® Safe-T-Pro
Uno, Assure® Lance
and Assure® Lance Plus



Side button
MPHSTSIDE21

Compares to
Unistik® 3



Pressure activated
MPHPRESS23

Compares to Microtainer®,
ReadyLance®, and
Mediance® Plus



Push button
MPHSFTY23

Compares to
Acti-Lance®



Pressure activated
MPHST28

Compares to SurgiLance®



Heel incision device
MPHNEWBORN2

Compares to Tenderfoot®
and Quikheel®

Safety Lancets

Easy for your staff, gentle for your patients

All of our lancets are designed for safe and easy punctures

Safety engineered

- The needle retracts automatically and locks, helping prevent accidental finger sticks
- Cannot be reloaded or reused, minimizing cross-contamination risk

Easy sampling

- Different gauges ensure the proper amount of blood is drawn
- Choose the style you prefer: push button, pressure activated or side button

Patient comfort

- Spring action and controlled penetration depth provide a quick, gentle puncture

General Purpose Lancet Items

Item No.	Description	Pkg.
MPH100030	General Purpose Lancets, 30G	2,500/cs
MPH100030Z	General Purpose Lancets, 30G	100/bx
MPHLANCINGPN2	Adjustable Lancing Device, Nonsafety	1/ea

Push Button

Item No.	Description	Pkg.
MPHSFTY28	28 gauge, 1.6 mm	3,600/cs
MPHSFTY28Z	28 gauge, 1.6 mm	200/bx
MPHSFTY281	28 gauge, 1.6 mm	2,000/cs
MPHSFTY281Z	28 gauge, 1.6 mm	100/bx
MPHSFTY285	28 gauge, 1.6 mm	1,200/cs
MPHSFTY285Z	28 gauge, 1.6 mm	50/bx
MPHSFTY23	23 gauge, 1.8 mm	2,000/cs
MPHSFTY23Z	23 gauge, 1.8 mm	200/bx
MPHSFTY21	21 gauge, 1.8 mm	3,600/cs
MPHSFTY21Z	21 gauge, 1.8 mm	200/bx
MPHSFTYBL	1.5 mm blade, 1.6 mm depth	3600/cs
MPHSFTYBLZ	1.5 mm blade, 1.6 mm depth	200/bx

Pressure Activated

Item No.	Description	Pkg.
MPHST28	28 gauge, 1.8 mm	2,000/cs
MPHST28Z	28 gauge, 1.8 mm	200/bx
MPHST21	21 gauge, 2.2 mm	2,000/cs
MPHST21Z	21 gauge, 2.2 mm	200/bx
MPHST21B	21 gauge, 1.8 mm	2,000/cs
MPHST21BZ	21 gauge, 1.8 mm	200/bx



Push button: Place the lancet firmly on chosen site and push the top button to activate



Side button: Place the lancet firmly on chosen site and push the side button to activate



Pressure activated (MPHST series): Push the lancet firmly against chosen site to activate



Pressure activated (MPHPRESS series): Push the lancet firmly against chosen site to activate

Side Button

Item No.	Description	Pkg.
MPHSTSIDE28	28 gauge, 1.8 mm	2,000/cs
MPHSTSIDE28Z	28 gauge, 1.8 mm	200/bx
MPHSTSIDE23	23 gauge, 1.8 mm	2,000/cs
MPHSTSIDE23Z	23 gauge, 1.8 mm	200/bx
MPHSTSIDE21	21 gauge, 2.4 mm	2,000/cs
MPHSTSIDE21Z	21 gauge, 2.4 mm	200/bx

Pressure Activated

- Expiration date printed on each individual lancet

Item No.	Description	Pkg.
MHPRESS30	30 gauge, 1.5 mm	2,000/cs
MHPRESS30Z	30 gauge, 1.5 mm	100/bx
MHPRESS28	28 gauge, 1.8 mm	2,000/cs
MHPRESS28Z	28 gauge, 1.8 mm	100/bx
MHPRESS26	26 gauge, 1.8 mm	2,000/cs
MHPRESS26Z	26 gauge, 1.8 mm	100/bx
MHPRESS23	23 gauge, 2.2 mm	2,000/cs
MHPRESS23Z	23 gauge, 2.2 mm	100/bx
MHPRESS21	21 gauge, 2.2 mm	2,000/cs
MHPRESS21Z	21 gauge, 2.2 mm	100/bx

Infant Heel Incision Devices

Handled with care

Designed with the precision you demand


Our easy-to-use trigger activates the blade quickly, safely and precisely making sure that each incision is made at the depth and width required for consistent blood flow


- High-quality stainless steel blade
- Small, precise incisions reduce bruising
- Color-coded sizing for easy identification
- Disposable, one time use

The defined contact area allows for accurate blade placement




Our integrated safety lock protects against accidental finger sticks and eliminates any loose pieces of plastic

Micro-preemie	Micro-Preemie	Item No.	Description	Pkg.
		MPHMPREEMIE2	Light Blue, 0.65 mm depth, 1.4 mm length	200/cs
		MPHMPREEMIE2Z	Light Blue, 0.65 mm depth, 1.4 mm length	50/bx

Preemie	Preemie	Item No.	Description	Pkg.
		MPHPREEMIE2	White, 0.85 mm depth, 1.75 mm length	200/cs
		MPHPREEMIE2Z	White, 0.85 mm depth, 1.75 mm length	50/bx

Newborn	Newborn	Item No.	Description	Pkg.
		MPHNEWBORN2	Pink, 1.0 mm depth, 2.5 mm length	200/cs
		MPHNEWBORN2Z	Pink, 1.0 mm depth, 2.5 mm length	50/bx

Toddler	Toddler	Item No.	Description	Pkg.
		MPHTODDLER2	Purple, 2.0 mm depth, 3.0 mm length	200/cs
		MPHTODDLER2Z	Purple, 2.0 mm depth, 3.0 mm length	50/bx

Increase blood flow

To ensure consistent blood flow, prep the baby's heel with an Infant Heel warmer to help minimize potential bruising, infection or cartilage damage caused by frequent heel sticks from blood tests



- Specially designed gel pack conforms gently to infant's heel
- Lift-and-stick tape for easy application and removal
- Activates instantly

Item No.	Description	Pkg.
MDS138007	Infant Heel Warmer, Gel Pack	100/cs

Total support

All your blood sampling needs in one place



NON21420



NON25650



MDS090735



OAT2585



MDS705153

Woven Gauze Sponges

- Folded edges prevent unraveling
- Both sides are packed in coated paper envelopes for easy opening

Item No.	Description	Pkg.
NON21420	2" X 2", 8-Ply, Sterile 2s	50 pk/bx, 30 bx/cs, 1,500 pks/cs
NON21424	4" X 4", 12-Ply, Sterile 2s	25 pk/bx, 24 bx/cs, 600 pks/cs

CURAD® Flex-Fabric™ Adhesive Bandage

- Woven bandage stretches and conforms to wound
- Fabric provides airflow to wound to help healing

Item No.	Description	Pkg.
NON25650	.75" X 3", Comfort Cloth	200/bx, 6 bx/cs

Alcohol Prep Pad

- Removes oils and residue from skin
- Soft, absorbent, non-woven pad is saturated with 70% isopropyl alcohol
- Intended for use on body only

Item No.	Description	Pkg.
MDS090735	Alcohol Prep Pad	200/bx, 3,000/cs

Restore® Sense™ Nitrile Exam Gloves

- Coated with layer of maxOat+®
- The National Eczema Association has awarded the Seal of Acceptance to Restore gloves

Item No.	Description	Pkg.
OAT2583	Restore Nitrile Glove, X-Small	2,500/cs
OAT2584	Restore Nitrile Glove, Small	2,500/cs
OAT2585	Restore Nitrile Glove, Medium	2,500/cs
OAT2586	Restore Nitrile Glove, Large	2,500/cs
OAT2587	Restore Nitrile Glove, X-Large	2,300/cs

Patient Room Sharps Containers

- For use in patient room, on med carts and in high-traffic areas
- Curved lid design facilitates one-handed disposal

Item No.	Description	Pkg.
MDS705153	5-Quart Sharps Container	20/cs
MDS705154	5-Quart Sharps Container	20/cs

Learn more. Contact your Medline Representative, call 1-800 MEDLINE (633-5463) or visit medline.com



We make healthcare run better™

Medline Industries, LP | 1-800-MEDLINE (633-5463) | www.medline.com

© 2025 Medline Industries, LP. All rights reserved. Medline is a registered trademark of Medline Industries, LP. All product names, regardless of whether or not they appear with a trademark symbol, are the property of their respective companies unless otherwise noted. MKT19W6740939 / e20909 / 94